

ENDOSCOPIC LINEAR CUTTER STAPLERS AND RELOADS

Endoscopic Linear Cutter Staplers are used in open and endoscopic surgeries for mechanical sutures, incisions and anastomoses. This version of Staplers is suitable for working with articulating and straight reloads with equal or different levels of staples. Endoscopic Linear Cutter Staplers are supplied in 3 different versions depending on suture length of 30, 45, 60 mm that suitable for tissues of various thicknesses: from vascular to super thick.

- Working part rotates 360° and allow to position reload at the required angle.
- Reload articulation 45° to each side.
- Available curved tip reload to provide visibility and maneuverability around target tissues and vessels.
- Anti-reflective coating of the shaft.
- Shaft diameter 12 mm which makes them suitable to work with trocars.
- 6 rows of compact ranged staples provide good hemostasia result.
- Knife is built in the reload.

Articulating Reloads with Equal Staple Height Size range

Tissue type/ reload color	Suture length, mm	Rows number	Open staple height, mm	Closed staple height, mm
VASCULAR	30	6	2.5	1
STANDARD	30	6	3.5	1.5
THICKENED	30	6	4.8	2
VASCULAR	45	6	2.5	1
STANDARD	45	6	3.5	1.5
SEMI-THICKENED	45	6	4.1	1.75
THICKENED	45	6	4.8	2
VASCULAR	60	6	2.5	1
STANDARD	60	6	3.5	1.5
SEMI-THICKENED	60	6	4.1	1.75
THICKENED	60	6	4.8	2



Find out more at www.invumed.com

INVUMED Innovations LTD. Cyprus. 23 John Kennedy Avenue «Globe House» - 1075 Nicosia
 invumed.com / info@invumed.com / +357 26220591
 © 2022 INVUMED Innovations LTD. All rights reserved. 27.12.22

Articulating Reloads with Different Staple Height Size range

Tissue type/ reload color	Suture length, mm	Rows number	Open staple height, mm	Closed staple height, mm
VASCULAR AND NORMAL	45	6	2.0/2.5/3.0	0.88 - 1.5
NORMAL AND THICKENED	45	6	3.0/3.5/4.0	1.5 - 2.25
SUPER-THICK	45	6	4.0/4.5/5.0	2.25 - 3.0
VASCULAR AND NORMAL	50	6	2.0/2.5/3.0	0.88 - 1.5
NORMAL AND THICKENED	50	6	3.0/3.5/4.0	1.5 - 2.25
SUPER-THICK	50	6	4.0/4.5/5.0	2.25 - 3.0

Rostral Articulating Reloads with Different Staple Height Size range

Tissue type/ reload color	Suture length, mm	Rows number	Open staple height, mm	Closed staple height, mm
VASCULAR AND NORMAL	45	6	2.0/2.5/3.0	0.88 - 1.5
NORMAL AND THICKENED	45	6	3.0/3.5/4.0	1.5 - 2.25
VASCULAR AND NORMAL	50	6	2.0/2.5/3.0	0.88 - 1.5
NORMAL AND THICKENED	50	6	3.0/3.5/4.0	1.5 - 2.25



Find out more at www.invumed.com

INVUMED Innovations LTD. Cyprus. 23 John Kennedy Avenue «Globe House» - 1075 Nicosia
 invumed.com / info@invumed.com / +357 26220591
 © 2022 INVUMED Innovations LTD. All rights reserved. 27.12.22

